

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Toxoplasma gondii IgM μ -capture**
 Prod.No.: **RE57111**
 Lot: **TOXM-244**
 Expiry: **2027-02-28**
 IFU: **Rev. 01**


Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0460-244	2025-08	2027-02-28
CONJS	1	KM0460-244	2025-08	2027-02-28
CONTROL -	1	NM0460-244	2025-08	2027-02-28
CUT OFF	1	CM0460-244	2025-08	2027-02-28
CONTROL +	1	PM0460-244	2025-08	2027-02-28
WASH BUF 20x	1	W-1117	2025-07	2027-07-31
SOLN STOP	1	S-1048	2025-05	2027-05-31
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-493	2025-05	2027-05-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.013	< Cut-off
CUT OFF	0.488	0.150-1.300
CONTROL +	1.266	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	25.84	14.21 - 37.47

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-09-05

Approved by: 
 Quality Control

Approved on: 2025-09-05



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany