



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Parvovirus B19 IgG**
 Prod.No.: **RE57281**
 Lot: **PARG-203**
 Expiry: **2027-04-30**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0370-203	2025-10	2027-04-30
CONJ	1	KG0370-203	2025-10	2027-04-30
CONTROL -	1	NG0370-203	2025-10	2027-04-30
CUT OFF	1	CG0370-203	2025-10	2027-04-30
CONTROL +	1	PG0370-203	2025-10	2027-04-30
WASH BUF 20x	1	W-1123	2025-10	2027-10-31
SOLN STOP	1	S-1054	2025-08	2027-08-31
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-495	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
CONTROL -	0.017	< 0.200 < Cut-off
CUT OFF	0.720	0.150-1.300
CONTROL +	1.803	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	29.99	16.49 - 43.49


The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-11-05

Approved by: 
 Quality Control

Approved on: 2025-11-05

 Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany

 IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany