

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Parvovirus B19 IgG**  
 Prod.No.: **RE57281**  
 Lot: **PARG-204**  
 Expiry: **2027-05-31**  
 IFU: **Rev. 01**


Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0370-204	2025-11	2027-05-31
<b>CONJ</b>	1	KG0370-204	2025-11	2027-05-31
<b>CONTROL -</b>	1	NG0370-204	2025-11	2027-05-31
<b>CUT OFF</b>	1	CG0370-204	2025-11	2027-05-31
<b>CONTROL +</b>	1	PG0370-204	2025-11	2027-05-31
<b>WASH BUF 20x</b>	1	W-1127	2025-11	2027-11-30
<b>SOLN STOP</b>	1	S-1062	2025-11	2027-11-30
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL</b>	1	D-496	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.012	< 0.200 < Cut-off
<b>CUT OFF</b>	0.563	0.150-1.300
<b>CONTROL +</b>	1.492	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	25.31	13.92 - 36.70


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-12-16

Approved by:   
 Quality Control

Approved on: 2025-12-16

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany