

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Parvovirus B19 IgM**

Prod.No.: **RE57291**

Lot: **PARM-202**

Expiry: **2027-03-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MM0370-202	2025-09	2027-03-31
<b>CONJ</b>	1	KM0370-202	2025-09	2027-03-31
<b>CONTROL -</b>	1	NM0370-202	2025-09	2027-03-31
<b>CUT OFF</b>	1	CM0370-202	2025-09	2027-03-31
<b>CONTROL +</b>	1	PM0370-202	2025-09	2027-03-31
<b>WASH BUF 20x</b>	1	W-1120	2025-08	2027-08-31
<b>SOLN STOP</b>	1	S-1051	2025-08	2027-08-31
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL M</b>	1	DM-347	2025-07	2027-07-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.009	< 0.200 < Cut-off
<b>CUT OFF</b>	0.552	0.150-1.300
<b>CONTROL +</b>	1.074	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	20.33	11.18 - 29.47

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-10-09

Approved by:   
Quality Control

Approved on: 2025-10-09



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany