


Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Varicella-Zoster Virus (VZV) IgG**
 Prod.No.: **RE57371**
 Lot: **VZVG-191**
 Expiry: **2027-04-30**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0490-191	2025-10	2027-04-30
CONJ	1	KG0490-191	2025-10	2027-04-30
CONTROL -	1	NG0490-191	2025-10	2027-04-30
CUT OFF	1	CG0490-191	2025-10	2027-04-30
CONTROL +	1	PG0490-191	2025-10	2027-04-30
WASH BUF 20x	1	W-1124	2025-10	2027-10-31
SOLN STOP	1	S-1055	2025-09	2027-09-30
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-495	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
CONTROL -	0.009	< 0.200 < Cut-off
CUT OFF	0.578	0.150-1.300
CONTROL +	1.155	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	21.18	11.65 - 30.70


The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-10-31

Approved by: 
 Quality Control

Approved on: 2025-10-31

 Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany

 IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany