



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Herpes simplex Virus 1 (HSV 1) IgM**
 Prod.No.: **RE57911**
 Lot: **HSV1M-144**
 Expiry: **2026-12-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0500-144	2025-06	2026-12-31
CONJ	1	KM0500-144	2025-06	2026-12-31
CONTROL -	1	NM0500-144	2025-06	2026-12-31
CUT OFF	1	CM0500-144	2025-06	2026-12-31
CONTROL +	1	PM0500-144	2025-06	2026-12-31
WASH BUF 20x	1	W-1113	2025-04	2027-04-30
SOLN STOP	1	S-1044	2025-03	2027-03-31
SUB TMB	1	T-078	2025-04	2028-04-28
DIL M	1	DM-345	2025-05	2027-05-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.018	< 0.200 < Cut-off
CUT OFF	0.477	0.150-1.300
CONTROL +	0.738	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	17.12	>10.0 - 24.82

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-07-10

Approved by: 
 Quality Control

Approved on: 2025-07-10



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany