

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Herpes simplex Virus 2 (HSV 2) IgG**  
 Prod.No.: **RE57921**  
 Lot: **HSV2G-144**  
 Expiry: **2027-06-30**  
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0540-144	2025-12	2027-06-30
<b>CONJ</b>	1	KG0540-144	2025-12	2027-06-30
<b>CONTROL -</b>	1	NG0540-144	2025-12	2027-06-30
<b>CUT OFF</b>	1	CG0540-144	2025-12	2027-06-30
<b>CONTROL +</b>	1	PG0540-144	2025-12	2027-06-30
<b>WASH BUF 20x</b>	1	W-1127	2025-11	2027-11-30
<b>SOLN STOP</b>	1	S-1062	2025-11	2027-11-30
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL</b>	1	D-496	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.015	< 0.200 < Cut-off
<b>CUT OFF</b>	0.454	0.150-1.300
<b>CONTROL +</b>	1.081	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	24.42	13.43 - 35.41


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-12-15

Approved by:   
 Quality Control

Approved on: 2025-12-15

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany