



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Herpes simplex Virus 2 (HSV 2) IgM**
 Prod.No.: **RE57931**
 Lot: **HSV2M-145**
 Expiry: **2027-08-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0540-145	2026-02	2027-08-31
CONJ	1	KM0540-145	2026-02	2027-08-31
CONTROL -	1	NM0540-145	2026-02	2027-08-31
CUT OFF	1	CM0540-145	2026-02	2027-08-31
CONTROL +	1	PM0540-145	2026-02	2027-08-31
WASH BUF 20x	1	W-1136	2026-02	2028-02-29
SOLN STOP	1	S-1068	2026-01	2028-01-31
SUB TMB	1	T-079	2025-12	2028-12-10
DIL M	1	DM-351	2026-01	2028-01-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.045	< 0.200 < Cut-off
CUT OFF	0.475	0.150-1.300
CONTROL +	0.923	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	19.77	10.87 - 28.66

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-03-17

Approved by: 
 Quality Control

Approved on: 2026-03-17



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany