



**GOLD  
STANDARD  
DIAGNOSTICS**



# Certificate of Analysis (CoA)

Product: **Herpes simplex Virus 1+2 (HSV) IgG**  
 Prod.No.: **RE57941**  
 Lot: **HSVG-121**  
 Expiry: **2027-09-30**  
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0250-121	2026-03	2027-09-30
<b>CONJ</b>	1	KG0250-121	2026-03	2027-09-30
<b>CONTROL -</b>	1	NG0250-121	2026-03	2027-09-30
<b>CUT OFF</b>	1	CG0250-121	2026-03	2027-09-30
<b>CONTROL +</b>	1	PG0250-121	2026-03	2027-09-30
<b>WASH BUF 20x</b>	1	W-1138	2026-03	2028-03-31
<b>SOLN STOP</b>	1	S-1070	2026-02	2028-02-29
<b>SUB TMB</b>	1	T-079	2025-12	2028-12-10
<b>DIL</b>	1	D-500	2026-02	2028-02-29

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.015	< 0.200 < Cut-off
<b>CUT OFF</b>	0.597	0.150-1.300
<b>CONTROL +</b>	1.057	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	19.64	10.80 - 28.48

The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-04-02

Approved by:   
 Quality Control

Approved on: 2026-04-02



Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany



IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany