



**GOLD  
STANDARD  
DIAGNOSTICS**



# Certificate of Analysis (CoA)

Product: **Entamoeba histolytica IgG**

Prod.No.: **RE58631**

Lot: **ENTG-088**

Expiry: **2027-02-28**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0140-088	2025-08	2027-02-28
<b>CONJ</b>	1	KG0140-088	2025-08	2027-02-28
<b>CONTROL -</b>	1	NG0140-088	2025-08	2027-02-28
<b>CUT OFF</b>	1	CG0140-088	2025-08	2027-02-28
<b>CONTROL +</b>	1	PG0140-088	2025-08	2027-02-28
<b>WASH BUF 20x</b>	1	W-1118	2025-08	2027-08-31
<b>SOLN STOP</b>	1	S-1049	2025-07	2027-07-31
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL</b>	1	D-493	2025-05	2027-05-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.014	< 0.200 < Cut-off
<b>CUT OFF</b>	0.483	0.150-1.300
<b>CONTROL +</b>	1.128	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	22.92	12.60 - 33.23


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-09-10

Approved by:   
 Quality Control

Approved on: 2025-09-10

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany