



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Leishmania infantum IgG**
 Prod.No.: **RE58641**
 Lot: **LEIG-098N**
 Expiry: **2026-08-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0310-098	2025-02	2026-08-31
CONJ	1	KG0310-098	2025-02	2026-08-31
CONTROL -	1	NG0310-098	2025-02	2026-08-31
CUT OFF	1	CG0310-098	2025-02	2026-08-31
CONTROL +	1	PG0310-098	2025-02	2026-08-31
WASH BUF 20x	1	W-1106	2025-01	2027-01-31
SOLN STOP	1	S-1037	2024-11	2026-11-30
SUB TMB	1	T-077	2022-09	2027-03-31
DIL	1	D-490	2024-12	2026-12-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.013	< 0.200 < Cut-off
CUT OFF	0.688	0.150-1.300
CONTROL +	1.482	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	24.96	13.73 - 36.20

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-03-11

Approved by: ACr
 Quality Control

Approved on: 2025-03-11



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany