

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Toxocara canis IgG**  
 Prod.No.: **RE58721**  
 Lot: **TOCG-103**  
 Expiry: **2027-07-31**  
 IFU: **Rev. 02**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0450-103	2026-01	2027-07-31
<b>CONJ</b>	1	KG0450-103	2026-01	2027-07-31
<b>CONTROL -</b>	1	NG0450-103	2026-01	2027-07-31
<b>CUT OFF</b>	1	CG0450-103	2026-01	2027-07-31
<b>CONTROL +</b>	1	PG0450-103	2026-01	2027-07-31
<b>WASH BUF 20x</b>	1	W-1131	2026-01	2028-01-31
<b>SOLN STOP</b>	1	S-1065	2025-12	2027-12-31
<b>SUB TMB</b>	1	T-079	2025-12	2028-12-10
<b>DIL</b>	1	D-498	2025-12	2027-12-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.004	< 0.200 < Cut-off
<b>CUT OFF</b>	0.426	0.150-1.300
<b>CONTROL +</b>	0.852	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	20.23	11.12 - 29.33


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-01-26

Approved by:   
 Quality Control

Approved on: 2026-01-26

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany