

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Taenia solium IgG**
 Prod.No.: **RE58731**
 Lot: **TAEG-084**
 Expiry: **2026-12-31**
 IFU: **Rev. 01**

| Component | Number | Lot | Prod. Date | Expiry Date |
|---------------------|--------|------------|------------|-------------|
| MTP | 1 | MG0420-084 | 2025-06 | 2026-12-31 |
| CONJ | 1 | KG0420-084 | 2025-06 | 2026-12-31 |
| CONTROL - | 1 | NG0420-084 | 2025-06 | 2026-12-31 |
| CUT OFF | 1 | CG0420-084 | 2025-06 | 2026-12-31 |
| CONTROL + | 1 | PG0420-084 | 2025-06 | 2026-12-31 |
| WASH BUF 20x | 1 | W-1113 | 2025-04 | 2027-04-30 |
| SOLN STOP | 1 | S-1044 | 2025-03 | 2027-03-31 |
| SUB TMB | 1 | T-078 | 2025-04 | 2028-04-28 |
| DIL | 1 | D-492 | 2025-03 | 2027-03-31 |

| Component | OD Value | OD Validation Criteria |
|------------------|----------|------------------------|
| CONTROL - | 0.017 | < 0.200 < Cut-off |
| CUT OFF | 0.543 | 0.150-1.300 |
| CONTROL + | 1.008 | > Cut-off |

| Component | Target Value [NTU] | Target Range [NTU] |
|------------------|--------------------|--------------------|
| CONTROL + | 18.90 | 10.40 - 27.41 |

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-07-17

Approved by: 
 Quality Control

Approved on: 2025-07-17



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany