

**GOLD  
STANDARD  
DIAGNOSTICS**



# Certificate of Analysis (CoA)

Product: **Epstein-Barr Virus (VCA) IgM**  
 Prod.No.: **RE58761**  
 Lot: **EBVM-215N**  
 Expiry: **2026-09-30**  
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MM0150-215	2025-03	2026-09-30
<b>CONJ</b>	1	KM0150-215	2025-03	2026-09-30
<b>CONTROL -</b>	1	NM0150-215	2025-03	2026-09-30
<b>CUT OFF</b>	1	CM0150-215	2025-03	2026-09-30
<b>CONTROL +</b>	1	PM0150-215	2025-03	2026-09-30
<b>WASH BUF 20x</b>	1	W-1108	2025-02	2027-02-28
<b>SOLN STOP</b>	1	S-1039	2025-01	2027-01-31
<b>SUB TMB</b>	1	T-077	2022-09	2027-03-31
<b>DIL M</b>	1	DM-343	2024-12	2026-12-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.010	< 0.200 < Cut-off
<b>CUT OFF</b>	0.450	0.150-1.300
<b>CONTROL +</b>	0.809	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	19.05	10.48 - 27.62


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-04-04

Approved by:   
 Quality Control

Approved on: 2025-04-04

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany