



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Mycoplasma pneumoniae IgM**
 Prod.No.: **RE58781**
 Lot: **MYCM-206**
 Expiry: **2027-10-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0350-206	2026-04	2027-10-31
CONJ	1	KM0350-206	2026-04	2027-10-31
CONTROL -	1	NM0350-206	2026-04	2027-10-31
CUT OFF	1	CM0350-206	2026-04	2027-10-31
CONTROL +	1	PM0350-206	2026-04	2027-10-31
WASH BUF 20x	1	W-1139	2026-03	2028-03-31
SOLN STOP	1	S-1071	2026-03	2028-03-31
SUB TMB	1	T-080	2026-02	2029-02-23
DIL M	1	DM-352	2026-03	2028-02-29

Component	OD Value	OD Validation Criteria
CONTROL -	0.097	< 0.200 < Cut-off
CUT OFF	0.494	0.150-1.300
CONTROL +	0.912	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	19.51	10.73 - 28.28

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-05-05

Approved by: Hoezalp
 Quality Control

Approved on: 2026-05-05



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany