



**GOLD
STANDARD
DIAGNOSTICS**



Certificate of Analysis (CoA)

Product: **Chikungunya Virus IgM μ -capture**
 Prod.No.: **RE58841**
 Lot: **CHIM-123**
 Expiry: **2027-05-31**
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0590-123	2025-11	2027-05-31
CONJ STRE	1	KM0590-123	2025-11	2027-05-31
AG LYO	6	AGM0590-123	2025-11	2027-05-31
AB SOLN	1	ABM0590-123	2025-11	2027-05-31
CONTROL -	1	NM0590-123	2025-11	2027-05-31
CUT OFF	1	CM0590-123	2025-11	2027-05-31
CONTROL +	1	PM0590-123	2025-11	2027-05-31
WASH BUF 20x	1	W-1128	2025-11	2027-11-30
SOLN STOP	1	S-1059	2025-10	2027-10-31
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-496	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
CONTROL -	0.244	< Cut-off
CUT OFF	0.610	0.150-1.300
CONTROL +	2.253	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	29.02	15.96 - 42.08

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-12-04



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany

Approved by: 
 Quality Control

Approved on: 2025-12-04



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany