

Quality Control Certificate



Fecal Helicobacter pylori-Antigen ELISA

REF RE58891


LOT AB641

2028-08-31
YYYY-MM-DD

Final Release Results:

	Calibrator concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.067	-	3%
CAL B	2.8	0.149	-	7%
CAL C	8.3	0.284	-	14%
CAL D	19.4	0.487	-	24%
CAL E	46.1	0.859	-	42%
CAL F	152.6	2.034	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	5.6	1.9	5.6	9.3
CONTROL 2	27.6	11.8	31.6	51.3



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-09-03

Accepted

This document has been produced electronically and is valid without a signature.