

Quality Control Certificate



Fecal Helicobacter pylori-Antigen ELISA

REF RE58891


LOT W979

2026-06-30
YYYY-MM-DD

Final Release Results:

	Calibrator concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.035	-	2%
CAL B	0.5	0.117	-	6%
CAL C	2.2	0.368	-	19%
CAL D	9.9	0.820	-	41%
CAL E	35.7	1.487	-	75%
CAL F	128.0	1.981	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	1.2	0.5	1.0	1.4
CONTROL 2	14.4	13.0	21.7	30.3



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2024-09-11

Accepted

This document has been produced electronically and is valid without a signature.