

Quality Control Certificate



Fecal Helicobacter pylori-Antigen ELISA

REF RE58891


LOT Y814

2027-02-28
YYYY-MM-DD

Final Release Results:

	Calibrator concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.082	-	3%
CAL B	2.6	0.213	-	8%
CAL C	7.2	0.448	-	17%
CAL D	24.9	0.768	-	29%
CAL E	49.4	1.239	-	47%
CAL F	161.1	2.652	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	3.8	2.5	4.6	6.7
CONTROL 2	32.2	18.7	31.2	43.6



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-12-17

Accepted

This document has been produced electronically and is valid without a signature.