

# Quality Control Certificate



## Fecal Helicobacter pylori-Antigen ELISA

REF RE58891


LOT Y992

2027-02-28  
YYYY-MM-DD

### Final Release Results:

	Calibrator concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.082	-	3%
CAL B	2.6	0.220	-	8%
CAL C	7.2	0.429	-	15%
CAL D	24.9	0.753	-	26%
CAL E	49.4	1.315	-	45%
CAL F	161.1	2.925	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	5.5	2.5	4.6	6.7
CONTROL 2	40.9	18.7	31.2	43.6



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-08-20

**Accepted**

This document has been produced electronically and is valid without a signature.