

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Leptospira IgG**

Prod.No.: **RE58931**

Lot: **LEPG-086N**

Expiry: **2026-07-31**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0660-086	2025-01	2026-07-31
<b>CONJ</b>	1	KG0660-086	2025-01	2026-07-31
<b>CONTROL -</b>	1	NG0660-086	2025-01	2026-07-31
<b>CUT OFF</b>	1	CG0660-086	2025-01	2026-07-31
<b>CONTROL +</b>	1	PG0660-086	2025-01	2026-07-31
<b>WASH BUF 20x</b>	1	W-1103	2024-10	2026-10-31
<b>SOLN STOP</b>	1	S-1034	2024-10	2026-10-31
<b>SUB TMB</b>	1	T-077	2022-09	2027-03-31
<b>DIL</b>	1	D-489	2024-10	2026-10-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.034	< 0.200 < Cut-off
<b>CUT OFF</b>	0.438	0.150-1.300
<b>CONTROL +</b>	0.753	> Cut-off

Final Report issued: 2025-02-10

Approved by: *HOe*  
Quality Control

Approved on: 2025-02-10



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany