

Quality Control Certificate



Serotonin ELISA


REF RE59121**LOT** ESE163**2027-07-31**

YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	manual	
			Range OD	OD/ODmax
CAL A	0	2.168	> 1.000	100%
CAL B	0.08	1.888	-	87%
CAL C	0.24	1.490	-	69%
CAL D	0.73	0.959	-	44%
CAL E	2.2	0.463	-	21%
CAL F	6.6	0.238	-	11%
CAL G	19.8	0.142	< 0.500	7%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	88	49	81	113
CONTROL 2	324	192	274	357



Kit Composition:

REF		LOT	
KESE651	LYO CONTROL 1	K875	2028-07
KESE661	LYO CONTROL 2	K876	2028-07
KEZZ811	MTP	K861	2029-06
KESE521	ANTISERUM	K872	2027-07
KESE531	BIOTIN	K871	2027-08
KESE541	CONC ENZCONJ	K877	2028-07
KESE611	CAL A-G	K870	2027-07

REF		LOT	
KESE761	ACYLREAG	K873	2027-07
KEZZ681	CONC ASSAYBUF	K833	2028-06
KEZZ701	CONC WASHBUF	K842	2028-06
KEZZ961	PNPP SUBS	K180-3	2027-09
KEZZ951	PNPP STOP	K874	2029-07

Valid edition of the instructions for use: V2025-01

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-08-10

Accepted

This document has been produced electronically and is valid without a signature.

**IBL International GmbH**
Flughafenstr. 52A, 22335 Hamburg, GermanyTel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl