

Quality Control Certificate



5-HIAA (5-Hydroxy-3-indole Acetic Acid) ELISA

REF RE59131


LOT EIN117

2027-09-30
YYYY-MM-DD

Final Release Results:

	Concentration	Measured value	Range	OD/ODmax
	mg/L	OD	OD	
CAL A	0.00	2.281	> 1.300	100%
CAL B	0.2	1.916	-	84%
CAL C	0.55	1.770	-	78%
CAL D	1.8	1.440	-	63%
CAL E	7.5	1.047	-	46%
CAL F	20	0.744	-	33%
CAL G	55	0.466	-	20%

	Concentration found	Acceptable range		
	mg/L	lower limit	target	upper limit
	mg/L	mg/L	mg/L	mg/L
CONTROL 1	4.3	2.7	3.9	5.1
CONTROL 2	27.8	20.2	28.1	35.9



Kit Composition:

REF		LOT	
KEIN651	CONTROL 1	K864	2028-01
KEIN661	CONTROL 2	K865	2028-01
KEZZ811	MTP	K828	2029-06
KEIN521	ANTISERUM	K218-A	2027-10
KEIN531	BIOTIN	K829	2027-10
KEIN541	CONC ENZCONJ	K863	2028-06
KEIN611	CAL A-G	K834	2027-10

REF		LOT	
KEZZ681	CONC ASSAYBUF	L950-D	2028-06
KEIN771	DILREAG	K832	2027-10
KEIN751	METHYLREAG	K830	2027-10
KEZZ701	CONC WASHBUF	K842	2028-06
KEZZ961	PNPP SUBS	K180-3	2027-09
KEZZ951	PNPP STOP	L737	2027-11
KEIN761	HCl	K831	2029-06

Valid edition of the instructions for use: V2014-11

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: S. Geweke

Validation date: 2026-08-07

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl