

Quality Control Certificate




Normetanephine ELISA

REF RE59171	LOT ENM202	2027-11-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration µg/L	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.537	> 1.000	100%
CAL B	77	2.123	-	84%
CAL C	192	1.677	-	66%
CAL D	480	1.062	-	42%
CAL E	1200	0.591	-	23%
CAL F	3000	0.309	-	12%
CAL G	7500	0.175	< 0.500	7%



	Concentration found µg/L	Acceptable range		
		lower limit µg/L	target µg/L	upper limit µg/L
CONTROL 1	224	133	208	284
CONTROL 2	642	427	667	907

Kit Composition:

REF		LOT	⌚
KEZZ621	CONTROL 1	K509	2027-12
KEZZ631	CONTROL 2	K510	2027-12
KEZZ811	MTP	K514	2028-12
KEZZ661	HYDRTUB	K491	2030-11
KEZZ671	HCl	L949	2028-04
KEZZ651	INDICATORBUF	K493	2028-11
KEZZ641	ACYLREAG	K507	2027-12
KENM541	CONC ENZCONJ	K527	2027-12

REF		LOT	⌚
KEZZ611	CAL A-G	K508	2027-12
KENM533	BIOTIN	K505	2027-12
KENM521	ANTISERUM	K506	2027-12
KEZZ681	CONC ASSAYBUF	L950-B	2027-12
KLZZ879	CONC WASHBUF	K494	2027-12
KEZZ961	PNPP SUBS	K492	2027-11
KEZZ951	PNPP STOP	K189	2028-05

Valid edition of the instructions for use: 2024-04

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-01-26

Accepted

This document has been produced electronically and is valid without a signature.