

Quality Control Certificate




Metanephrine ELISA

REF RE59181	LOT EMN207	2027-11-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration µg/L	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.780	> 1.000	100%
CAL B	26	2.503	-	90%
CAL C	64	2.199	-	79%
CAL D	160	1.756	-	63%
CAL E	400	1.268	-	46%
CAL F	1000	0.735	-	26%
CAL G	2500	0.407	< 0.500	15%

	Concentration found µg/L	Acceptable range		
		lower limit µg/L	target µg/L	upper limit µg/L
CONTROL 1	115.1	77.6	121.2	164.8
CONTROL 2	583	363	568	772



Kit Composition:

REF		LOT	⌚
KEZZ621	CONTROL 1	K577	2028-01
KEZZ631	CONTROL 2	K578	2028-01
KEZZ811	MTP	K634	2029-02
KEZZ661	HYDRTUB	K537	2030-12
KEZZ671	HCL	K580	2029-01
KEZZ651	INDICATORBUF	K493	2028-11
KEZZ641	ACYLREAG	K579	2028-01
KEMN541	CONC ENZCONJ	K641	2028-02

REF		LOT	⌚
KEZZ611	CAL A-G	K576	2028-01
KEMN533	BIOTIN	K639	2028-02
KEMN521	ANTISERUM	K638	2028-02
KEZZ681	CONC ASSAYBUF	L950-C	2028-01
KLZZ879	CONC WASHBUF	K582	2028-01
KEZZ961	PNPP SUBS	K492	2027-11
KEZZ951	PNPP STOP	L737	2027-11

Valid edition of the instructions for use: 2024-04

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-03-25

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl