

# Quality Control Certificate




## Histamine ELISA (Plasma)

<b>REF</b> RE59221	<b>LOT</b> EHP201	<b>2026-07-31</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

### Final Release Results:

	Concentration ng/mL	Measured value OD	Range	OD/ODmax
CAL A	0	2.292	> 1.000	100%
CAL B	0.20	2.024	-	88%
CAL C	1.2	1.663	-	73%
CAL D	8.0	1.063	-	46%
CAL E	35	0.554	-	24%
CAL F	120	0.232	< 0.500	10%



	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	3.3	1.3	2.6	3.9
CONTROL 2	19.6	10.6	17.7	24.8

### Kit Composition:

REF		LOT	
KEZZ811	MTP	L940	2028-03
KEHP521	ANTISERUM	L929	2026-07
KEHP541	CONC ENZCONJ	L953	2026-08
KEHP611	CAL P A-F	L931	2026-10
KEHP651	CONTROL P 1	L932	2026-10
KEHP661	CONTROL P 2	L933	2026-10
KEHP751	ACYLREAG	L934	2026-07

REF		LOT	
KEHP761	CONC ASSAYBUF	L930	2027-03
KEZZ701	CONC WASHBUF	L939	2027-03
KEZZ651	INDICATORBUF	L792	2027-12
KEZZ991	TMB SUBS	L919	2028-06
KEZZ921	TMB STOP	L296	2027-04

Valid edition of the instructions for use: 2024-03

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-05-07

**Accepted**

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)