

# Quality Control Certificate



## Histamine ELISA (Urine)

REF RE59221


LOT EHP201

2026-07-31  
YYYY-MM-DD

### Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.204	> 1.000	100%
CAL B	2.7	1.985	-	90%
CAL C	8.1	1.706	-	77%
CAL D	24.3	1.237	-	56%
CAL E	73	0.614	-	28%
CAL F	219	0.340	< 0.500	15%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	11.2	7.1	10.8	14.4
CONTROL 2	39.0	32.5	41.7	50.8



### Kit Composition:

REF		LOT	
KEZZ811	MTP	L940	2028-03
KEHP521	ANTISERUM	L929	2026-07
KEHP541	CONC ENZCONJ	L953	2026-08
KEHP621A	CAL U/C A	N483-6	2026-07
KEHP621	CAL U/C B-F	N484E-5D	2026-07
KEHP671	CONTROL U/C 1	L216A-2A	2026-07
KEHP681	CONTROL U/C 2	L217A-2A	2026-07

REF		LOT	
KEHP751	ACYLREAG	L934	2026-07
KEHP761	CONC ASSAYBUF	L930	2027-03
KEZZ701	CONC WASHBUF	L939	2027-03
KEZZ651	INDICATORBUF	L792	2027-12
KEZZ991	TMB SUBS	L919	2028-06
KEZZ921	TMB STOP	L296	2027-04

Valid edition of the instructions for use: 2024-03

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-05-07

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)