

Quality Control Certificate



Histamine ELISA (Plasma)

REF RE59221


LOT EHP202

2027-03-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range	OD/ODmax
CAL A	0	1.998	> 1.000	100%
CAL B	0.20	1.634	-	82%
CAL C	1.2	1.295	-	65%
CAL D	8.0	0.831	-	42%
CAL E	35	0.410	-	21%
CAL F	120	0.175	< 0.500	9%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	1.8	1.2	2.5	3.7
CONTROL 2	11.7	8.0	13.3	18.6



Kit Composition:

REF		LOT	
KEZZ811	MTP	K442	2028-10
KEHP521	ANTISERUM	K443	2027-04
KEHP541	CONC ENZCONJ	K444	2027-04
KEHP611	CAL P A-F	K445	2027-06
KEHP651	CONTROL P 1	K446A	2027-06
KEHP661	CONTROL P 2	K447A	2027-06
KEHP751	ACYLREAG	K452	2027-03

REF		LOT	
KEHP761	CONC ASSAYBUF	K453	2027-11
KEZZ701	CONC WASHBUF	K280-2	2027-07
KEZZ651	INDICATORBUF	K464	2028-11
KEZZ991	TMB SUBS	K132-7	2028-09
KEZZ921	TMB STOP	K441	2028-11

Valid edition of the instructions for use: 2025-07

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-02-09

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl