

Quality Control Certificate



Histamine ELISA (Urine)

REF RE59221


LOT EHP202

2027-03-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.105	> 1.000	100%
CAL B	2.7	1.793	-	85%
CAL C	8.1	1.513	-	72%
CAL D	24.3	1.034	-	49%
CAL E	73	0.561	-	27%
CAL F	219	0.268	< 0.500	13%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	8.4	5.9	8.9	11.9
CONTROL 2	38.2	29.1	37.3	45.5



Kit Composition:

REF		LOT	
KEZZ811	MTP	K442	2028-10
KEHP521	ANTISERUM	K443	2027-04
KEHP541	CONC ENZCONJ	K444	2027-04
KEHP621A	CAL U/C A	N483-7	2027-03
KEHP621	CAL U/C B-F	K478	2027-03
KEHP671	CONTROL U/C 1	K474	2027-03
KEHP681	CONTROL U/C 2	K475	2027-03

REF		LOT	
KEHP751	ACYLREAG	K452	2027-03
KEHP761	CONC ASSAYBUF	K453	2027-11
KEZZ701	CONC WASHBUF	K280-2	2027-07
KEZZ651	INDICATORBUF	K464	2028-11
KEZZ991	TMB SUBS	K132-7	2028-09
KEZZ921	TMB STOP	K441	2028-11

Valid edition of the instructions for use: 2025-07

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-02-09

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl