

Quality Control Certificate




CatCombi (Adrenalin/Noradrenalin) ELISA (Noradrenalin)

REF RE59242	LOT ECC186	2026-06-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.108	-	5%
CAL B	5.0	0.196	-	9%
CAL C	15	0.379	-	17%
CAL D	50	0.894	-	39%
CAL E	150	1.675	-	73%
CAL F	500	2.279	-	100%

	Concentration found ng/mL	Acceptable range		
	ng/mL	lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	21.5	14.4	21.2	28.0
CONTROL 2	118	78	115	152



Kit Composition:

REF	LOT		REF	LOT	
KEZZ521	CONTROL 1	K212 2027-06	KEZZ575	LYO COMT	K204 2026-06
KEZZ531	CONTROL 2	K213 2027-06	KEZZ581	COENZ	K203 2026-06
KEZZ801	MTP	K209 2028-06	KEZZ591	ENZBUF	K201 2028-06
KEZZ411	ANTISERUM NAD	K206 2026-07	KEZZ422	COMT ADD	K202 2026-07
KEZZ821	CONC ENZCONJ	K214 2027-07	KEZZ601/9	RELEASEBUF	K200 2026-07
KEZZ511	CAL A-F	K211 2027-06	KLZZ879	CONC WASHBUF	K215 2027-06
KEZZ541	ACYLREAG	K208 2026-07	KEZZ971	PNPP SUBS	L747-4 2026-08
KEZZ555	EXTRPLATE	K210 2027-06	KEZZ951	PNPP STOP	K216 2028-06
KEZZ565	EXTRBUF	K199 2026-07			

Valid edition of the instructions for use: 2014-12

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T.Schmidt

Validation date: 2025-09-01

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl