

Quality Control Certificate




Noradrenalin (Norepinephrine) ELISA

| | | |
|--------------------|-------------------|---------------------------------|
| REF RE59261 | LOT END184 | 2026-06-30 YYYY-MM-DD |
|--------------------|-------------------|---------------------------------|

Final Release Results:

| | Concentration ng/mL | Measured value OD | Range OD | OD/ODmax |
|-------|------------------------|-------------------------|-------------|----------|
| CAL A | 0 | 0.108 | - | 5% |
| CAL B | 5.0 | 0.196 | - | 9% |
| CAL C | 15 | 0.379 | - | 17% |
| CAL D | 50 | 0.894 | - | 39% |
| CAL E | 150 | 1.675 | - | 73% |
| CAL F | 500 | 2.279 | - | 100% |

| | Concentration found ng/mL | Acceptable range | | |
|-----------|---------------------------------|----------------------|-----------------|----------------------|
| | | lower limit ng/mL | target ng/mL | upper limit ng/mL |
| CONTROL 1 | 21.5 | 14.4 | 21.2 | 28.0 |
| CONTROL 2 | 118 | 78 | 115 | 152 |



Kit Composition:

| REF | LOT | 2026-06 | REF | LOT | 2026-06 | | |
|---------|---------------|---------|---------|-----------|--------------|--------|---------|
| KEZZ521 | CONTROL 1 | K212 | 2027-06 | KEZZ575 | LYO COMT | K204 | 2026-06 |
| KEZZ531 | CONTROL 2 | K213 | 2027-06 | KEZZ581 | COENZ | K203 | 2026-06 |
| KEZZ801 | MTP | K209 | 2028-06 | KEZZ591 | ENZBUF | K201 | 2028-06 |
| KEZZ411 | ANTISERUM NAD | K206 | 2026-07 | KEZZ422 | COMT ADD | K202 | 2026-07 |
| KEZZ821 | CONC ENZCONJ | K214 | 2027-07 | KEZZ601/9 | RELEASEBUF | K200 | 2026-07 |
| KEZZ511 | CAL A-F | K211 | 2027-06 | KLZZ879 | CONC WASHBUF | K215 | 2027-06 |
| KEZZ541 | ACYLREAG | K208 | 2026-07 | KEZZ971 | PNPP SUBS | L747-4 | 2026-08 |
| KEZZ555 | EXTRPLATE | K210 | 2027-06 | KEZZ951 | PNPP STOP | K216 | 2028-06 |
| KEZZ565 | EXTRBUF | K199 | 2026-07 | | | | |

Valid edition of the instructions for use: 2014-12

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T.Schmidt

Validation date: 2025-09-01

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl