

Quality Control Certificate



Neopterin ELISA

REF RE59321


LOT ENO241

2027-05-31
YYYY-MM-DD

Final Release Results:

	Concentration nmol/L	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.557	> 1.000	100%
CAL B	1.35	2.080	-	81%
CAL C	4	1.542	-	60%
CAL D	12	0.953	-	37%
CAL E	37	0.510	-	20%
CAL F	111	0.257	-	10%

	Concentration found nmol/L	Acceptable range		
		lower limit nmol/L	target nmol/L	upper limit nmol/L
CONTROL 1	5.6	3.4	5.6	7.9
CONTROL 2	19.1	11.8	18.2	24.5



Kit Composition:

REF		LOT	
KENO651	CONTROL 1	K127	2027-05
KENO661	CONTROL 2	K128	2027-05
KEZZ811	MTP	K124	2028-05
KENO521	ANTISERUM	K125	2027-05
KENO541	ENZCONJ	K129	2027-06
KENO611	CAL A-F	K126	2027-05
KENO751	ASSAYBUF	K130	2027-05
KEZZ751	CONC WASHBUF	K131	2027-06

REF		LOT	
KEZZ981	TMB SUBS	K132	2028-09
KENO781	TMB STOP	K133	2028-05

Valid edition of the instructions for use: 2024-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: M.Dries

Validation date: 2025-07-30

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl