

Quality Control Certificate




Neopterin ELISA

REF RE59321	LOT ENO242	2028-01-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration nmol/L	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.597	> 1.000	100%
CAL B	1.35	2.081	-	80%
CAL C	4	1.538	-	59%
CAL D	12	0.950	-	37%
CAL E	37	0.495	-	19%
CAL F	111	0.248	-	10%

	Concentration found nmol/L	Acceptable range		
		lower limit nmol/L	target nmol/L	upper limit nmol/L
CONTROL 1	5.9	3.6	5.9	8.3
CONTROL 2	21.9	13.8	21.3	28.7



Kit Composition:

REF		LOT	
KENO651	CONTROL 1	K499	2028-01
KENO661	CONTROL 2	K500	2028-01
KEZZ811	MTP	K495	2029-01
KENO521	ANTISERUM	K496	2028-01
KENO541	ENZCONJ	K497	2028-01
KENO611	CAL A-F	K498	2028-01
KENO751	ASSAYBUF	K501	2028-01
KEZZ751	CONC WASHBUF	K502	2028-01

REF		LOT	
KEZZ981	TMB SUBS	K503	2029-04
KENO781	TMB STOP	K504	2028-12

Valid edition of the instructions for use: 2024-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-02-19

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl