

Quality Control Certificate



FinTest™ IgG4 Screen 90 ELISA

REF**RE59499-292220****LOT****IK90-145****2028-07-21**

YYYY-MM-DD

Final Release Results:

	Concentration U/mL	Result OD	Acceptable range OD
CAL A	0.35	0.105	≥ 0.040
CAL B	0.7	0.192	≥ 0.075
CAL C	3.5	0.322	≥ 0.125
CAL D	17.5	0.788	≥ 0.300
CAL E	50	1.628	≥ 0.575
CAL F	100	2.083	≥ 0.750

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Jörg Spangenberg

Accepted

Validation date: 2026-09-03

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl