

Quality Control Certificate




Amyloid-beta (1-40) CSF ELISA

REF RE59651	LOT EAY129	2026-12-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	0.009	< 0.080	0%
CAL B	118	0.092		5%
CAL C	235	0.212		11%
CAL D	470	0.517		26%
CAL E	940	0.992		50%
CAL F	1880	1.987	-	100%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	757	490	700	910
CONTROL 2	179	124	177	230



Kit Composition:

REF		LOT	
KEAY511	MTP	L985	2027-03
KEAY611	CAL A	L986	2027-03
KEAY621	CAL B-F	L987	2027-03
KEAY651	CONTROL 1	L988	2027-03
KEAY661	CONTROL 2	L989	2027-03

REF		LOT	
KEAY541	CONC ENZCONJ	L990	2027-03
KEAY571	ASSAYBUF	L873	2026-12
KEZZ451	SAMPLEDIL	L872	2026-12
KEZZ441	CONC WASHBUF	L938	2027-02
KEZZ991	TMB SUBS	L541-6	2027-10
KEZZ921	TMB STOP	L861	2028-01

Valid edition of the instructions for use: 2019-05

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-04-17

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl