

# Quality Control Certificate




## Progesterone Luminescence Immunoassay

<b>REF</b> RE62021 & RE62029	<b>LOT</b> LPR159	<b>2026-12-31</b> YYYY-MM-DD
------------------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration pg/mL	Measured value RLU	Range RLU	RLU/RLU <sub>max</sub>
CAL A	0	1063365	-	100%
CAL B	10	960165	-	90%
CAL C	25	842528	-	79%
CAL D	50	661743	-	62%
CAL E	100	488768	-	46%
CAL F	300	243448	-	23%
CAL G	1000	110413	-	10%



	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	46.4	29.7	45.8	61.8
CONTROL 2	186	149	186	223

**Kit Composition:**

REF		LOT	
KLPR651/9	CONTROL 1	K163	2026-12
KLPR661/9	CONTROL 2	K164	2026-12
KLZZ511	MTP	K155	2027-11
KLPR541/9	CONC ENZCONJ	K187	2027-01
KLPR611/9	CAL A	K161	2026-12
KLPR621/9	CAL B-G	K162	2026-12
KLPR521/9	ANTISERUM	K174	2026-12

REF		LOT	
KLZZ681/9	ASSAYBUF	K156	2027-05
KLZZ601/9	LUMINREAG AP	K115-2	2027-03
KLZZ879	CONC WASHBUF	K301	2027-07

Valid edition of the instructions for use: V2018-07

**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: M.Dries

**Accepted**

Validation date: 2025-08-27

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)