

Quality Control Certificate




DHEA Luminescence Immunoassay

REF RE62051 & RE62059	LOT LDH147	2026-09-30 YYYY-MM-DD
------------------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	3724810	-	100%
CAL B	12.3	2581630	-	69%
CAL C	37	2091060	-	56%
CAL D	111	1577250	-	42%
CAL E	333	1114988	-	30%
CAL F	1000	697565	-	19%
CAL G	3000	425528	-	11%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	83	37	91	146
CONTROL 2	455	356	475	594



Kit Composition:

REF	LOT	Icon
KLDH651/9	CONTROL 1	L909
KLDH661/9	CONTROL 2	L910
KLDH511	MTP	M792
KLDH541/9	CONC ENZCONJ	L907
KLDH611/9	CAL A	L908
KLDH621/9	CAL B-G	L908
KLDH521/9	ANTISERUM	L906

REF	LOT	Icon
KLZZ681/9	ASSAYBUF	L721
KLDH561/9	LUMINREAG AP	K115
KLZZ879	CONC WASHBUF	L885

Valid edition of the instructions for use: V2015_01

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: **Dr. Spangenberg**

Validation date: **2025-05-23**

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl