

Quality Control Certificate




Cortisol Saliva Luminescence Immunoassay

REF RE62111 & RE62119	LOT LCL204	2026-02-28 YYYY-MM-DD
------------------------------	-------------------	---------------------------------

Final Release Results:

	Concentration µg/dL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	30476355	-	100%
CAL B	0.015	22161403	-	73%
CAL C	0.03	17782108	-	58%
CAL D	0.10	9197910	-	30%
CAL E	0.40	2915055	-	10%
CAL F	1.00	1162213	-	4%
CAL G	3.20	310425	-	1%

	Concentration found µg/dL	Acceptable range		
		lower limit µg/dL	target µg/dL	upper limit µg/dL
CONTROL 1	0.080	0.051	0.079	0.107
CONTROL 2	0.72	0.51	0.73	0.95



Kit Composition:

REF		LOT	
KLCL651/9	CONTROL 1	L488	2026-08
KLCL661/9	CONTROL 2	L489	2026-08
KLCL511	MTP	L403	2026-08
KLCL541/9	ENZCONJ	L494	2026-08
KLCL611/9	CAL A	L486	2026-08
KLCL621/9	CAL B-G	L487	2026-08
KLZZ709	CONC WASHBUF	L490	2026-08

REF		LOT	
KLZZ641/9	LUMINREAG 1	L243-2	2026-02
KLZZ651/9	LUMINREAG 2	L244-2	2026-02

Valid edition of the instructions for use: 2022-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2024-08-06

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl