

Quality Control Certificate




17beta-Estradiol Luminescence Immunoassay

REF RE62141	LOT LEL139	2027-04-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	1703485	-	100%
CAL B	0.9	1589040	-	93%
CAL C	2	1488160	-	87%
CAL D	4	1191433	-	70%
CAL E	8	972440	-	57%
CAL F	16	462078	-	27%
CAL G	64	206158	-	12%



	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	3.1	2.0	3.1	4.1
CONTROL 2	13.3	11.0	13.7	16.5

Kit Composition:

REF		LOT	
KLEL651	CONTROL 1	K344	2027-04
KLEL661	CONTROL 2	K345	2027-04
KLEL511	MTP	L918	2027-08
KLEL541	CONC ENZCONJ	K417	2027-05
KLEL611	CAL A	K343	2027-04
KLEL621	CAL B-G	K343	2027-04

REF		LOT	
KLEL521	ANTISERUM	K342	2027-04
KLZZ681	ASSAYBUF	K341	2027-08
KLZZ601	LUMINREAG AP	K416	2027-08
KLZZ879	CONC WASHBUF	K388	2027-10

Valid edition of the instructions for use: 2023-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-11-26

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl