

Quality Control Certificate




17beta-Estradiol Luminescence Immunoassay

REF RE62141	LOT LEL140	2027-04-30 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	2374435	-	100%
CAL B	0.9	2221778	-	94%
CAL C	2	2023858	-	85%
CAL D	4	1769405	-	75%
CAL E	8	1342118	-	57%
CAL F	16	637210	-	27%
CAL G	64	289583	-	12%



	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	3.4	2.1	3.3	4.4
CONTROL 2	13.0	10.2	12.8	15.4

Kit Composition:

REF		LOT	
KLEL651	CONTROL 1	K457	2027-05
KLEL661	CONTROL 2	K458	2027-05
KLEL511	MTP	K340	2028-02
KLEL541	CONC ENZCONJ	K418	2027-05
KLEL611	CAL A	K456	2027-05
KLEL621	CAL B-G	K456	2027-05

REF		LOT	
KLEL521	ANTISERUM	K342	2027-04
KLZZ681	ASSAYBUF	K341	2027-08
KLZZ601	LUMINREAG AP	K416	2027-08
KLZZ879	CONC WASHBUF	K388	2027-10

Valid edition of the instructions for use: 2023-09

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Accepted

Validation date: 2025-11-27

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl