

Quality Control Certificate




IRT Neonatal Luminescence Immunoassay

REF RE68019	LOT LIT200	2027-01-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	42725	-	2%
CAL B	25	142905	-	8%
CAL C	48	304340	-	16%
CAL D	95	521260	-	28%
CAL E	185	1105335	-	59%
CAL F	370	1870670	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	67	54	78	101
CONTROL 2	119	88	125	163



Kit Composition:

REF		LOT	⌚	REF		LOT	⌚
KLIT655	CONTROL 1	256	2027-01	KLIT755	CONC ASSAYBUF	L884	2027-02
KLIT665	CONTROL 2	257	2027-01	KLZZ635	LUMINREAG	L917	2027-02
KLIT511	MTP	L883	2027-02	KLZZ879	CONC WASHBUF	L885	2027-02
KLIT545	CONC ENZCONJ	L915	2027-02				
KLIT615	CAL A-F	255	2027-01				

Valid edition of the instructions for use: 2024-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: M.Dries

Accepted

Validation date: 2025-03-03

This document has been produced electronically and is valid without a signature.