

Quality Control Certificate




IRT Neonatal Luminescence Immunoassay

REF RE68019	LOT LIT201	2027-02-28 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	31710	-	2%
CAL B	25	99855	-	6%
CAL C	48	173690	-	11%
CAL D	95	467210	-	29%
CAL E	185	783460	-	49%
CAL F	370	1589690	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	80	56	80	104
CONTROL 2	138	99	141	184



Kit Composition:

REF	LOT		REF	LOT	
KLIT655	CONTROL 1	262	2027-02	KLIT755	CONC ASSAYBUF
KLIT665	CONTROL 2	263	2027-02	KLZZ635	LUMINREAG
KLIT511	MTP	L883	2027-02	KLZZ879	CONC WASHBUF
KLIT545	CONC ENZCONJ	L941	2027-03		
KLIT615	CAL A-F	261	2027-02		

Valid edition of the instructions for use: 2024-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-04-09

This document has been produced electronically and is valid without a signature.