

Quality Control Certificate




IRT Neonatal Luminescence Immunoassay

REF RE68019	LOT LIT204	2027-10-31 YYYY-MM-DD
--------------------	-------------------	---------------------------------

Final Release Results:

	Concentration ng/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	43770	-	2%
CAL B	20	166045	-	7%
CAL C	45	255545	-	11%
CAL D	90	745705	-	32%
CAL E	180	1360555	-	58%
CAL F	400	2360520	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target ng/mL	upper limit ng/mL
CONTROL 1	71	53	76	99
CONTROL 2	116	89	127	165



Kit Composition:

REF		LOT	⌚	REF		LOT	⌚
KLIT655	CONTROL 1	285	2027-10	KLIT755	CONC ASSAYBUF	K469	2027-11
KLIT665	CONTROL 2	286	2027-10	KLZZ635	LUMINREAG	K470	2027-10
KLIT511	MTP	K460	2027-11	KLZZ879	CONC WASHBUF	K388	2027-10
KLIT545	CONC ENZCONJ	K473	2027-11				
KLIT615	CAL A-F	284	2027-10				

Valid edition of the instructions for use: 2024-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Accepted

Validation date: 2025-12-17

This document has been produced electronically and is valid without a signature.