


Quality Control Certificate



Cardiolipin IgG ELISA

REF RE75011

LOT 25070115

 2026-06-30
YYYY-MM-DD

Final Release Results:

	Conc. GPL-U/mL	Measured value OD	Acceptable range		
			lower limit OD	target value OD	upper limit OD
CAL A	0	0.068			< 0.100
CAL F	100	1.972	> 1.000		
CONTROL -	-	0.105			≤ 0.300
CONTROL +	-	0.910	≥ 0.500		≤ 1.600

	Validation criteria (manual)			
	Measured value GPL-U/mL	Acceptable range		
		lower limit GPL-U/mL	target value GPL-U/mL	upper limit GPL-U/mL
CONTROL -	1.2			≤ 6.0
CONTROL +	30	19	30	41

Factor for qualitative assay evaluation	0.48
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-07-28

This document has been produced electronically and is valid without a signature.