


Quality Control Certificate



Jo-1-Ab ELISA

REF RE75271

LOT 25020147

 2026-05-31
YYYY-MM-DD

Final Release Results:

	Conc. U/mL	Measured value OD	Acceptable range		
			lower limit OD	target value OD	upper limit OD
CAL A	0	0.051			< 0.100
CAL F	100	2.036	> 1.000		
CONTROL -	-	0.118			≤ 0.300
CONTROL +	-	1.117	≥ 0.500		≤ 1.600

	Validation criteria (manual)			
	Measured value U/mL	Acceptable range		
		lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL -	0.5			≤ 2.0
CONTROL +	21	17	27	37

Factor for qualitative assay evaluation	0.28
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-03-13

This document has been produced electronically and is valid without a signature.