


# Quality Control Certificate



## Jo-1-Ab ELISA

REF RE75271

LOT 26030170

 2027-07-31  
YYYY-MM-DD

### Final Release Results:

	Conc. U/mL	Measured value OD	Acceptable range		
			lower limit OD	target value OD	upper limit OD
CAL A	0	0.055			< 0.100
CAL F	100	2.246	> 1.000		
CONTROL -	-	0.145			≤ 0.300
CONTROL +	-	1.428	≥ 0.500		≤ 1.600

	Validation criteria (manual)			
	Measured value U/mL	Acceptable range		
		lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL -	0.6			≤ 2.0
CONTROL +	25	17	27	37

Factor for qualitative assay evaluation	0.28
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

# Accepted

Validation date: 2026-05-21

This document has been produced electronically and is valid without a signature.