


Quality Control Certificate



Modified Gliadin Peptide IgA ELISA

REF RE75741

LOT 25020133

 2026-07-31
YYYY-MM-DD

Final Release Results:

	Conc.	Measured value	Acceptable range		
	U/mL	OD	lower limit OD	target value OD	upper limit OD
Calibrator A	0	0.075			< 0.100
Calibrator F	100	2.501	> 1.000		
CONTROL -	-	0.151			≤ 0.300
CONTROL +	-	1.240	≥ 0.500		≤ 1.600

	Validation criteria (manual)			
	Measured value	Acceptable range		
	U/mL	lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL -	2.0			≤ 5.0
CONTROL +	39	23	36	49

Factor for qualitative assay evaluation	0.28
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2025-03-21

This document has been produced electronically and is valid without a signature.