

Quality Control Certificate



Modified Gliadin Peptide IgG ELISA

REF RE75751

LOT 25080231

 2027-02-28
YYYY-MM-DD

Final Release Results:

	Conc.	Measured value	Acceptable range		
	U/mL	OD	lower limit OD	target value OD	upper limit OD
Calibrator A	0	0.054	> 1.000		
Calibrator F	100	2.158			
CONTROL -	-	0.128	≥ 0.500		
CONTROL +	-	1.145			

	Validation criteria (manual)			
	Measured value	Acceptable range		
	U/mL	lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL -	1.6	23	36	≤ 5.0
CONTROL +	32			49

Factor for qualitative assay evaluation	0.31
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-10-07

Accepted

This document has been produced electronically and is valid without a signature.