


# Quality Control Certificate



## Modified Gliadin Peptide IgG ELISA

REF RE75751

LOT 26050067

 2027-10-31  
YYYY-MM-DD

### Final Release Results:

|              | Conc. | Measured value | Acceptable range  |                    |                   |
|--------------|-------|----------------|-------------------|--------------------|-------------------|
|              | U/mL  | OD             | lower limit<br>OD | target value<br>OD | upper limit<br>OD |
| Calibrator A | 0     | 0.046          |                   |                    | < 0.100           |
| Calibrator F | 100   | 2.472          | > 1.000           |                    |                   |
| CONTROL -    | -     | 0.157          |                   |                    | ≤ 0.300           |
| CONTROL +    | -     | 1.425          | ≥ 0.500           |                    | ≤ 1.600           |

|           | Validation criteria (manual) |                     |                      |                     |
|-----------|------------------------------|---------------------|----------------------|---------------------|
|           | Measured value<br>U/mL       | lower limit<br>U/mL | target value<br>U/mL | upper limit<br>U/mL |
| CONTROL - | 2.0                          |                     |                      | ≤ 5.0               |
| CONTROL + | 35                           | 23                  | 36                   | 49                  |

|   |      |
|---|------|
| Factor for qualitative assay evaluation | 0.31 |
|---|------|

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

# Accepted

Validation date: 2026-05-21

This document has been produced electronically and is valid without a signature.