

# Quality Control Certificate



## Phenylalanine (PKU) neonatal Screening Assay

<b>REF</b> RE80015 & RE80019	<b>LOT</b> NP231	<b>2022-02-28</b> YYYY-MM-DD
------------------------------	------------------	---------------------------------

### Final Release Results:

	Concentration mg/dL	Measured value OD	Range OD	OD/ODmax
CAL A	1.60	0.056	-	20%
CAL B	3.50	0.086	-	30%
CAL C	6.50	0.125	-	44%
CAL D	12.00	0.231	-	81%
CAL E	17.00	0.285	-	100%

	Concentration found mg/dL	Acceptable range		
		lower limit mg/dL	target value mg/dL	upper limit mg/dL
CONTROL 1	3.1	2.5	3.2	3.8
CONTROL 2	5.4	4.5	5.7	6.8

### Kit Composition:

<b>REF</b>		<b>LOT</b>	
KENP655	CONTROL 1	114	2022-02
KENP665	CONTROL 2	114	2022-02
KENP545	LYO ENZ	112	2022-08
KENP755/9	ENZDIL	991	2022-05
KENP615	CAL A-E	113	2022-02
KENP775/9	DYE	992	2022-05
KENP765/9	SUBS	990	2022-05

Valid edition of the instructions for use: 2019-02

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: S. Geweke

Validation date: 2021-08-26

**Accepted**

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)